



# Preventive maintenance

## Respironics V200 and Esprit Critical Care Ventilator service agreements

Preventive Maintenance Service Agreements support customers who desire to repair their own Respironics V200 or Esprit Ventilators but don't have time to do the preventive maintenance service. You can be confident that your ventilators will be maintained on schedule by product experts.

### Annual preventive maintenance service

- Installation of annual preventive maintenance parts kit
- Complete performance verification
- Electrical safety test
- Interior and exterior cleaning
- Latest base software revision (upon request)
- Visual inspection
- Up to \$100 (list price) in minor cosmetic repair

### 12,500-hour preventive maintenance service

- Installation of 12,500 hour preventive maintenance kit
- Complete performance verification
- Electrical safety test
- Interior and exterior cleaning
- Latest base software revision (upon request)
- Visual inspection
- Up to \$100 (list price) in minor cosmetic repairs

Items not covered by a Preventive Maintenance Service Agreement or Per-Call Preventive Maintenance Service include travel and accessories external to the ventilator such as, but not limited to carts, filters, oxygen sensors, oxygen hoses, patient circuit components, humidifiers, humidifier accessories, back-up batteries and external batteries

### Preventive-Maintenance-Only Service Agreements

Period	Part Number	Features
1 year	1041682	One on-site annual preventive maintenance service
2 years	1041681	Two on-site annual preventive maintenance services
1 year	1041685	One on-site 12,500-hour preventive maintenance service
2 years	1041684	Two on-site 12,500-hour preventive maintenance services*

\* If a second 12,500-hour preventive maintenance is not required due to low use, an annual preventive maintenance will be substituted.

### Per-Call Preventive Maintenance Services (non-agreement)

Period	Part Number	Features
1 year	1007501	Annual preventive maintenance service
2 years	1009540	12,500-hour preventive maintenance service
3 years	1001733	12,500-hour preventive maintenance kit (customers receive a return credit when the blower is returned to Philips)

**PHILIPS**



© 2009 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)  
fax: +31 40 27 64 887

MCI 4102338  
4522 962 51748 \* MAY 2009

Philips Healthcare  
Global Information Center  
P.O. Box 1286  
5602 BG Eindhoven  
The Netherlands